

# $^{41}\text{Ca}$ and Accelerator Mass Spectrometry to Monitor Calcium Metabolism in End Stage Renal Disease Patients

ROBERT L. FITZGERALD,<sup>1,2\*</sup> DARREN J. HILLEGONDS,<sup>3</sup> DOUGLAS W. BURTON,<sup>1</sup>  
TERRANCE L. GRIFFIN,<sup>1</sup> SCOTT MULLANEY,<sup>1,4</sup> JOHN S. VOGEL,<sup>3,5</sup> LEONARD J. DEFTOS,<sup>1,4</sup>  
and DAVID A. HEROLD<sup>1,2</sup>

**Background:** Monitoring bone resorption with measurements of bone density and biochemical markers is indirect. We hypothesized that bone resorption can be studied directly by serial measurements of the ratio  $^{41}\text{Ca}/\text{Ca}$  in serum after in vivo labeling of calcium pools with  $^{41}\text{Ca}$ . We report the preparation of an intravenous  $^{41}\text{Ca}$  dose suitable for humans, an analytical method for determining  $^{41}\text{Ca}/\text{Ca}$  isotope ratios in biological samples, and studies in human volunteers.

**Methods:**  $^{41}\text{Ca}$  was formulated and aliquoted into individual vials, and to the extent possible, the  $^{41}\text{Ca}$  doses were tested according to US Pharmacopeia (USP) guidelines. A 10 nCi dose of  $^{41}\text{Ca}$  was administered intravenously to 4 end stage renal disease (ESRD) patients on hemodialysis and 4 healthy control individuals. Distribution kinetics were determined over 168 days. Calcium was isolated with 3 precipitation steps and a cation-exchange column, and  $^{41}\text{Ca}/\text{Ca}$  ratios in serum were then measured by accelerator mass spectrometry.

**Results:** The dosing solution was chemically and radiologically pure, contained  $<0.1$  endotoxin unit/mL, and passed USP sterility tests. Quantification of  $^{41}\text{Ca}/\text{Ca}$  ratios was linear from  $6 \times 10^{-14}$  to  $9.1 \times 10^{-10}$ . The run-to-run imprecision (as CV) of the method was 4% at  $4.6 \times 10^{-11}$  and 6% at  $9.1 \times 10^{-10}$ . The area under the curve of  $^{41}\text{Ca}$  in the central compartment vs time was

significantly less for ESRD patients than for controls ( $P < 0.005$ ).

**Conclusions:** Isotope ratios spanning 5 orders of magnitude can be measured by accelerator mass spectrometry with excellent precision in the range observed in samples collected from patients who have received 10 nCi of  $^{41}\text{Ca}$ . The  $^{41}\text{Ca}$  at this dose caused no adverse effects in 8 volunteers. This is the first report of the use of  $^{41}\text{Ca}$  to monitor differences in bone turnover between healthy individuals and ESRD patients.

© 2005 American Association for Clinical Chemistry

Calcium exists naturally as 6 stable isotopes, with the most abundant isotope being  $^{40}\text{Ca}$ . Stable ( $^{42}\text{Ca}$ ,  $^{44}\text{Ca}$ , and  $^{46}\text{Ca}$ ) and short-lived radioactive isotopes ( $^{45}\text{Ca}$  and  $^{47}\text{Ca}$ ) of calcium have been used to study calcium metabolism for years (1–3). Limitations of the stable isotopes are that they are expensive to purify and that relatively large doses need to be given to distinguish the exogenously administered tracer from the naturally occurring signal. Limitations of the short-lived radioisotopes are that they deliver significant amounts of radiation and that their short half-lives prevent long-term tracing. The long-lived radioisotope,  $^{41}\text{Ca}$ , is an ideal candidate for biological tracer studies because the background signal of  $^{41}\text{Ca}$  is essentially zero as a result of the very low natural production of  $^{41}\text{Ca}$ . The advantage of combining a marker such as  $^{41}\text{Ca}$  (low noise) with a technique with attomole sensitivity, such as accelerator mass spectrometry (AMS),<sup>6</sup> is that excellent signal/noise ratios can be obtained from nanogram quantities of isotopes administered (4–12).

$^{41}\text{Ca}$  has a half-life of 104 000 years and decays by

<sup>1</sup> VA San Diego Healthcare System, San Diego, CA.

Departments of <sup>2</sup> Pathology and <sup>4</sup> Medicine, University of California, San Diego, CA.

<sup>3</sup> Lawrence Livermore National Laboratories, Livermore, CA.

<sup>5</sup> Department of Nutrition, University of California, Davis, CA.

\*Address correspondence to this author at: VAMC-113, 3350 La Jolla Village, San Diego, CA 92161. Fax 858-642-6304; e-mail rlfitzgerald@vapop.ucsd.edu.

Received February 14, 2005; accepted August 15, 2005.

Previously published online at DOI: 10.1373/clinchem.2005.049650

<sup>6</sup> Nonstandard abbreviations: AMS, accelerator mass spectrometry; ESRD, end stage renal disease; AP, alkaline phosphatase; USP, United States Pharmacopeia; and PTH, parathyroid hormone.

electron capture to the naturally occurring and stable isotope  $^{41}\text{K}$ , emitting x-rays and electrons of very low energy ( $<3.6$  keV) in the process (13, 14). For biological experiments,  $^{41}\text{Ca}$  is typically prepared by irradiating calcium samples with neutrons in a nuclear reactor. Other radioisotopes of calcium,  $^{45}\text{Ca}$  and  $^{47}\text{Ca}$ , are also formed during this process, but these have relatively short half-lives (163 and 4.5 days, respectively). After sufficient time is allowed for the  $^{45}\text{Ca}$  to decay, the only radioactive material remaining is  $^{41}\text{Ca}$ . The long lifetime of  $^{41}\text{Ca}$  prevents its routine analysis by decay counting, but AMS is particularly well suited for long-term monitoring of bone resorption by direct quantification of  $^{41}\text{Ca}$  concentrations. The long half-life of  $^{41}\text{Ca}$ , combined with the sensitivity of AMS, allows isotope ratios to be measured for a period of years, a time span that is not feasible with other isotopes of calcium. Potential clinical uses of  $^{41}\text{Ca}$  include monitoring of bone loss in patients with end stage renal disease (ESRD), diagnosis of osteoporosis, monitoring of the effectiveness of antiresorptive drugs, and early detection of bone metastasis.

Biochemical markers of bone disease have been in use for more than 30 years, but although beneficial in certain settings, all have drawbacks, primarily related to analytical and biological variability (15–19). One method for interpreting bone markers involves looking for the “least significant change”, which incorporates both the analytical and biological variability. The least significant change amounts for bone alkaline phosphatase (AP), osteocalcin, pyridinoline, deoxypyridinoline, and N-telopeptide are in the range of 15% to 40% (20). Bone mineral density, as measured by dual energy x-ray absorptiometry, provides a quantitative measure of the mineralization of bone but is relatively insensitive to change, generally taking 12 to 18 months before changes can be measured. A widely used diagnostic method for identifying bone disease is biopsy and histomorphometry; however, this procedure has morbidity and, in the absence of severe bone disease, is not clinically indicated (21). We report a novel technique for directly assessing the long-term turnover of the mineral phase of bone by measuring  $^{41}\text{Ca}/\text{Ca}$  ratios in serum after an intravenous dose of 10 nCi of  $^{41}\text{Ca}$ .

### Materials and Methods

In keeping with general standards used in isotope abundance work, we report  $^{41}\text{Ca}/\text{Ca}$ , where Ca reflects all isotopes of calcium.

#### PREPARATION OF THE DOSING MATERIAL

Initial purification and isotopic analysis of the  $^{41}\text{Ca}$  material used for the dosing solution have been described (22). We used an aliquot of the  $^{41}\text{Ca}$  master stock solution (described as “C1” in the original publication) as the starting material for the dosing solution (22). Radiologic purity is virtually guaranteed from the lack of significant elemental contaminants and the long decay period (the original dose material was purchased in 1984). However,

to ensure that we had no radioactive contaminants, we measured the emitted radiation with a gamma/beta survey meter and a high-purity germanium detector.

We received 2.060 g of a solution of  $\text{CaCO}_3$  in 4 mol/L  $\text{HNO}_3$  containing 0.7034 mg of calcium per gram of solution. By weight, the  $^{41}\text{Ca}$  abundance in the solution was 1.232% (22); consequently, this solution contained 18  $\mu\text{g}$  of  $^{41}\text{Ca}$  or 1500 nCi. The target for the dosing solution was 10 nCi of  $^{41}\text{Ca}$  in each vial as a 9-mL solution. In a sterile hood, using sterile technique, we added  $\sim 500$  mL of 0.9% sodium chloride injection (9 g/L) USP (Baxter Healthcare) and 3.75 mL of 10% USP  $\text{CaCl}_2$  (0.9 mol/L; American Reagent Labs) to a 2-L sterile vessel. The 2.060 g of  $^{41}\text{Ca}$  solution was quantitatively transferred to this vessel, and the pH was adjusted to 5–7 by addition of 7.5% sodium bicarbonate USP (0.9 mol/L; American Reagent Labs). The remainder of the solution was diluted to 1350 mL by weight. In the sterile hood, the dosing solution was filtered through a 0.2  $\mu\text{m}$  filter (Corning). With the solution still in the sterile hood, we volumetrically pipetted 9-mL portions into 10-mL clear vials (Wheaton Scientific). A gray straight plug (Wheaton Scientific) and aluminum seal (Wheaton Scientific) were placed loosely on the vial and held in place with a piece of autoclave indicating paper. The doses were autoclaved at 121  $^\circ\text{C}$  and 20 psi (140 kPa) for 30 min. In addition to the autoclave indicating tape, the autoclave process was monitored with a biological indicator (3M Attest) to determine that sterilization was successful. Immediately after autoclaving, the aluminum seals were crimped in the sterile hood. Representative samples were analyzed for endotoxins and sterility as described below. Throughout the process, plastic containers and pipetting devices were used to decrease the chance that trace metals would be introduced into the dosing solution. Pyrogenicity (endotoxin testing), sterility, bacterial, and fungistasis testing were then performed according to United States Pharmacopeia (USP) guidelines (23, 24). Elemental composition (purity) was checked by inductively coupled plasma/mass spectrometry using good manufacturing procedures (GMP). Inductively coupled plasma/mass spectrometry was chosen for purity testing because this method has excellent sensitivity ( $<1$   $\mu\text{g}/\text{L}$  for most elements), is broad spectrum (able to detect over 60 elements), and is the type of analysis best suited to identify potential impurities in the  $^{41}\text{Ca}$  dosing solution. There are no USP guidelines for acceptable limits of impurities for  $^{41}\text{Ca}$  dosing solutions; we therefore chose a conservative limit of  $<0.00001\%$  for any heavy metal contaminant [i.e., the maximum amount of heavy metals to which individuals would be exposed would be 0.00009 g (90  $\mu\text{g}$ )].

Radiologic purity of the  $^{41}\text{Ca}$  master stock solution (1.5  $\mu\text{Ci}$ ) was determined with a gamma/beta survey meter and with a high-purity germanium detector. To obtain an energy spectrum of the electron capture decay of  $^{41}\text{Ca}$  to  $^{41}\text{K}$ , the entire 1.5  $\mu\text{Ci}$  of  $^{41}\text{Ca}$  was placed in a high-purity germanium detector for 1 week.

## PATIENTS

The human subject protocol was approved by Institutional Review Boards at the University of California's Human Research Protections Program and Lawrence Livermore National Laboratory. ESRD patients were selected as the study group because 75%–100% of these patients have bone disease (25), and this study was designed as a proof-of-principle investigation. Our aim was to determine whether  $^{41}\text{Ca}/\text{Ca}$  ratios changed more rapidly in patients with known bone disease than in controls without known bone disease. To further enhance the differences between our control and study groups, we chose patients with high-turnover bone disease as evidenced by increased concentrations of parathyroid hormone (PTH). Because the absorption of ingested calcium in ESRD patients is known to vary widely, we chose to administer the dose intravenously to ensure that all controls and patients received the same initial dose of  $^{41}\text{Ca}$ .

Eight human volunteers participated in the study, 4 controls and 4 ESRD patients. All volunteers gave written informed consent. Vital signs and pain rating (0–10) of the volunteers were obtained before administration of 10 nCi of  $^{41}\text{Ca}$  via a heparin lock. To further ensure sterility, the entire 9 mL of solution to be administered was aspirated through a 0.22  $\mu\text{m}$  pore-size sterile filter into a 10-mL syringe. The sterile filter was removed, and the dose was administered by slow intravenous push through a heparin lock. The heparin lock was then flushed with 10 mL of normal saline. The heparin lock was removed, and vital signs and pain scores were obtained after administration. The ESRD patients received hemodialysis 3 times weekly with 3.5–4 h per treatment. Their dialysis was performed with high-flux membranes, and their calcium management followed K-DOQI guidelines (26).  $^{41}\text{Ca}$  was administered immediately after dialysis treatment.

## SAMPLE ANALYSIS

Samples were analyzed at the Center for AMS at Lawrence Livermore National Laboratory, which performs high-throughput, routine analyses of multiple nuclides, including  $^{41}\text{Ca}$  (4, 10). The automated AMS control system and streamlined sample preparation allow preparation and measurement of nuclides in more than 100 samples per day by a single operator.

Total AP, urea nitrogen, creatinine, potassium, sodium, chloride, carbon dioxide, phosphate, calcium, albumin, and glucose concentrations were measured with Beckman Synchron reagents on an LX20Pro analyzer. The PTH midmolecule was measured with a previously described assay (27, 28). 25-Hydroxyvitamin D, 1,25-dihydroxyvitamin D, PTH (intact; detects amino acids 7–84), N-telopeptide, and bone-specific AP were measured by ARUP Laboratories (Salt Lake City, UT).

Analysis of  $^{41}\text{Ca}$  in plasma by AMS was adapted from methods published by Lin et al. (29) and Freeman et al. (30). Elimination of the microwave digestion step simpli-

fied the analysis.  $^{41}\text{Ca}/\text{Ca}$  ratios were measured by AMS after calcium was isolated and purified using the steps described below.

We added 130  $\mu\text{L}$  of 30 g/L (0.75 mol/L) calcium as  $\text{CaCO}_3$  to 1 mL of plasma samples to ensure that there was sufficient calcium for AMS ( $\sim 4$  mg of Ca, or 0.1 mmol). Samples were then acid-digested for 30 min at room temperature with 2 mL of concentrated  $\text{HNO}_3$ . After digestion, 7 mL of deionized  $\text{H}_2\text{O}$  was added to the specimens; a yellow precipitant formed and was discarded after centrifugation (750g for 5 min), and the supernatant was retained. The calcium in the supernatant was then precipitated overnight with 1 mL of saturated ammonium oxalate and 3 mL of concentrated  $\text{NH}_4\text{OH}$ . Samples were centrifuged (750g for 10 min), the supernatant was discarded, and the resulting calcium oxalate was dissolved in 0.4 mL of 5 mol/L  $\text{HNO}_3$ . The samples were then diluted with deionized  $\text{H}_2\text{O}$  to an acid strength of 0.08 mol/L and chromatographed on a cation-exchange column as described below.

Cation-exchange columns (1.6 g of Bio-Rad AG 50W-X8, 200–400 mesh, hydrogen form) were prepared by washing with 5 mL of 5 mol/L  $\text{HNO}_3$ , 10 mL of deionized  $\text{H}_2\text{O}$ , and 5 mL of 0.08 mol/L  $\text{HNO}_3$ . The acid-digested plasma extracts were applied to the columns, which were then sequentially washed with two 4.5-mL portions of 0.08 mol/L  $\text{HNO}_3$ . The calcium was then eluted with two 4.25-mL portions of 5 mol/L  $\text{HNO}_3$ . The pH of the resulting solution was increased with 1.5 mL of deionized  $\text{H}_2\text{O}$  and 1.5 mL of concentrated ammonium hydroxide, and samples were allowed to cool for 30 min. After the samples were cool, 3 mL of concentrated HF was added, and  $\text{CaF}_2$  was allowed to precipitate overnight. (Note: concentrated HF is extremely lethal, and all users must be properly trained before using this acid.) Samples were centrifuged, and the resulting  $\text{CaF}_2$  was washed 3 times with 1 mL of deionized  $\text{H}_2\text{O}$ . After the  $\text{H}_2\text{O}$  was decanted, samples were placed in a 100  $^\circ\text{C}$  oven to dry for at least 12 h. Recovery of the added calcium was determined gravimetrically. Samples were then mixed with silver powder (2–4 parts  $\text{CaF}_2$  to silver by weight) to improve conductivity in the AMS source and loaded into aluminum targets for AMS analysis.

As outlined in Fig. 1, AMS  $^{41}\text{Ca}$  measurements were accomplished via (a) ionization of  $\text{CaF}_2$ /silver mixture, (b) selection of 97 or 98  $m/z$  ions (target ions are  $^{40}\text{CaF}_3^-$  and  $^{41}\text{CaF}_3^-$ , respectively), (c) acceleration of ions through  $9 \times 10^6$  V, (d) removal of the 3 fluorine atoms and 9 electrons from the calcium atom through charge exchange in a thin carbon foil at the high-voltage terminal to destroy any molecular interferences, (e) a second stage of acceleration to ground potential, (f) measurement of  $^{40}\text{Ca}^{8+}$  in an offset Faraday cup, and (g) ion identification and single-particle counting of  $^{41}\text{Ca}^{8+}$  in a multianode gas ionization detector. The details of the AMS analysis can be found elsewhere (4). Our analytical  $^{41}\text{Ca}/\text{Ca}$  background value was  $\sim 6 \times 10^{-14}$ . Measurement precision was monitored by




Fig. 1. Schematic of Lawrence Livermore accelerator mass spectrometer.

use of secondary standards similar in isotope ratio to unknowns and was <3%, assuming sufficient sample <sup>41</sup>Ca content to identify at least 1000 <sup>41</sup>Ca ions.

The method was validated at <sup>41</sup>Ca/Ca isotope ratios varying from blanks (~6 × 10<sup>-14</sup>) to concentrations in the range expected 2 h post dose (9.08 × 10<sup>-10</sup>). We added <sup>41</sup>Ca to these sera to achieve the target isotope ratios and processed them as described above. Within-run variability was determined at 6 isotope ratios (from blanks to 9.08 × 10<sup>-10</sup>) with n = 3 at each concentration. Run-to-run variability was determined at 3 isotope ratios (from blanks to 9.08 × 10<sup>-10</sup>) with each batch of samples. Run-to-run variability represented a minimum of 7 different batches over a 6-month period.

Results

Because the <sup>41</sup>Ca solution was given to human subjects, considerable effort was made to characterize it in terms of elemental composition, radioactivity, sterility, and pyrogenicity. The administered solution was free from heavy metal contaminants, and only 2 elements, sodium and calcium, were found at >1 part per million (Fig. 2); both had been added as part of the formulation process.

The radiologic purity of the <sup>41</sup>Ca master stock solution (1.5 μCi) was determined for the concentrated stock solution, rather than the diluted <sup>41</sup>Ca solution that was

given to volunteers, to ensure maximized sensitivity and radiologic purity. In the master stock solution, neither γ nor β radiation above background was detected (Table 1). In fact, according to US Nuclear Regulatory Commission guidelines, the measured results of the prepared dosing solution fell below the limit of what is considered radioactive (2 nCi/g). As can be seen in Fig. 3, purity testing gave the expected low-energy internal bremsstrahlung emissions from the electron capture decay of <sup>41</sup>Ca and no measurable β or γ activity, thus demonstrating that the <sup>41</sup>Ca stock solution was radiologically pure. The maximum energy of the sample occurred at ~421 keV, a result consistent with the electron capture decay of <sup>41</sup>Ca to <sup>41</sup>K (31).

Results of USP sterility tests showed no growth after 14 days, and USP endotoxin testing demonstrated <0.1 endotoxin units/mL. Taken together, these results demonstrate that the <sup>41</sup>Ca dose was radiologically and chemically pure and, importantly, also sterile and free of endotoxin.

The difference of results from expected values and the within-run imprecision of the method for <sup>41</sup>Ca/Ca ratios in serum samples with added <sup>41</sup>Ca over 4 orders of magnitude of isotope ratios expected when administering 10-nCi doses of <sup>41</sup>Ca to human volunteers are shown in Table 2. Run-to-run imprecision was in the 5% range with errors of 1%–3% (Table 3). As can also be seen in Table 3, with isotope ratios near background (~6 × 10<sup>-14</sup>), it is difficult to routinely obtain good accuracy and precision because the background amount of <sup>41</sup>Ca is so small that it is difficult to detect, an ideal situation for labeling studies. However, after a single 10-nCi dose, volunteers had <sup>41</sup>Ca/Ca ratios several orders of magnitude greater than background; therefore, the poor precision and accuracy at background isotope ratios was not relevant to our studies. The mean recovery of calcium through the purification process was 87%.

There were no postdose changes in vital signs or pain ratings, compared with predose values, for any of the 8 volunteers. The mean baseline clinical chemistry values




Fig. 2. Multielemental analysis of the administered <sup>41</sup>Ca solution.

White boxes indicate that elements were not analyzed; gray boxes represent elements that were present at <1 part per million; black boxes represent elements present at concentrations >1 part per million.



**Table 1. Result of gamma/beta measurement of stock  $^{41}\text{Ca}$  solution (1.5  $\mu\text{Ci}$ ) showing no detectable counts above blank with 2 separate gamma/beta counters.**

Instrument	Last calibrated	Calibration due	Mean (SD) counts	
			Blank	$^{41}\text{Ca}$ stock
E120 A	January 2003	April 2004	45 (15)	45 (15)
E120 B	September 2002	December 2003	30 (10)	30 (10)

for the controls and the ESRD patients are shown in Table 4. The patients were hyperparathyroid, and the increased concentrations of N-telopeptide, bone-specific AP, and total AP provide evidence of increased bone turnover.

The initial distribution and elimination kinetics of the 10-nCi  $^{41}\text{Ca}$  dose in the plasma from the controls (mean age, 54 years;  $n = 4$ ) and the ESRD patients [mean age, 59 years; mean (SD) time on dialysis, 2.0 (1.1) years;  $n = 4$ ] are shown in Fig. 4. As can be seen in Fig. 4, over the initial 28 days, the  $^{41}\text{Ca}/\text{Ca}$  isotope ratio in plasma decreased more rapidly in the ESRD patients than in the controls, giving a lower area under the curve for ESRD patients compared with the controls ( $P < 0.005$ ).

### Discussion

Any time humans are exposed to radiation, the principle of ALARA (as low as reasonably achievable) applies. We administered 10 nCi (120 ng) of  $^{41}\text{Ca}$  to 8 human subjects, 4 of whom were seriously ill. The biological risks associ-




Fig. 3. Acquired energy spectrum for 2.06 g (18.0  $\mu\text{g}$  of  $^{41}\text{Ca}$  = 1.5  $\mu\text{Ci}$ ) stock dose solution as measured for 1 week by use of a high-purity germanium detector.

No ionizing radiation (as  $\alpha$ ,  $\beta$ , or  $\gamma$  radiation) is produced; the low energy shown is from the electron capture decay of  $^{41}\text{Ca}$  to  $^{41}\text{K}$ .

**Table 2. Within-run precision ( $n = 3$ ) of  $^{41}\text{Ca}/\text{Ca}$  measurements in serum samples.**

Expected $^{41}\text{Ca}/\text{Ca}$	Measured $^{41}\text{Ca}/\text{Ca}^a$	CV, %	Error, %
$6.00 \times 10^{-14}$	$6.17 (0.30) \times 10^{-14}$	5	3
$4.55 \times 10^{-13}$	$5.01 (0.24) \times 10^{-14}$	5	10
$4.57 \times 10^{-11}$	$4.91 (0.17) \times 10^{-11}$	3	7
$9.14 \times 10^{-11}$	$9.61 (0.27) \times 10^{-11}$	3	5
$4.99 \times 10^{-10}$	$5.01 (0.04) \times 10^{-10}$	1	<1
$9.08 \times 10^{-10}$	$9.08 (0.05) \times 10^{-10}$	1	<1

<sup>a</sup> Mean (SD);  $n = 3$  at each concentration.

ated with radiation exposure at this low exposure are too small to quantify. The only other published report of administration of intravenous  $^{41}\text{Ca}$  to a human subject calculated that the overall radiation dose commitment was 0.06  $\mu\text{Sv}$  for the first year after a 10-nCi dose of  $^{41}\text{Ca}$  and that this value is 30 000 times smaller than the total dose from natural radiation over this period (32). In simple terms, the amount of radiation to which our volunteers were exposed is equivalent to that received from a 5-min commercial airplane ride (33,34). The amount of radioactivity that we administered was  $\sim 1000$ -fold less than that used in studies involving other calcium radioisotopes (35).

In addition to the perceived risk of radioactivity, other limitations of using  $^{41}\text{Ca}/\text{Ca}$  ratios to monitor bone health include the high cost of instrumentation and the extensive sample preparation required for analysis. The transition of AMS from the research laboratory to a clinical laboratory may not occur in the near future, but it is worth remembering that the first mass spectrometers, with very limited capabilities, began as large complicated instruments relegated to sophisticated research laboratories. Today, many clinical laboratories have several mass spectrometers with a variety of ionization and analyzer configurations. At present, there is sufficient AMS beam time to meet current needs. If the methods we are developing become essential for monitoring bone health, additional resources will need to be developed. Alternative techniques for measuring  $^{41}\text{Ca}$  include resonance ionization mass spectrometry (36) and atom trap trace analysis (37), but these methods have not yet matched the throughput or sensitivity of AMS.

Before  $^{41}\text{Ca}$  is used clinically, it must be demonstrated that measurement of the isotope ratio after a  $^{41}\text{Ca}$  dose has clinical relevance. Our study is the first to show that  $^{41}\text{Ca}/\text{Ca}$  ratios in a diseased group are different from

**Table 3. Interassay imprecision of  $^{41}\text{Ca}/\text{Ca}$  measurements.**

Expected $^{41}\text{Ca}/\text{Ca}$	Measured $^{41}\text{Ca}/\text{Ca}^a$	CV, %	Error, %	n
$6.00 \times 10^{-14}$	$5.59 (2.78) \times 10^{-14}$	50	7	22
$4.57 \times 10^{-11}$	$4.70 (0.19) \times 10^{-11}$	4	3	10
$9.13 \times 10^{-11}$	$9.20 (0.39) \times 10^{-11}$	4	1	18
$9.08 \times 10^{-10}$	$8.81 (0.52) \times 10^{-10}$	6	3	12

<sup>a</sup> Mean (SD).

**Table 4. Mean (SD) baseline laboratory values for controls and ESRD patients.**

Clinical test	Reference interval	Controls	ESRD patients	P
25-OH-vitamin D, ng/L	20–57	22 (8)	23 (15)	NS <sup>a</sup>
1,25-Vitamin D, ng/L	15–75	43 (18)	14 (5)	<0.05
PTH (intact), ng/L	15–75	40 (9)	800 (250)	<0.001
PTH (midmolecule), ng/L	5–70	13 (9)	1500 (750)	<0.01
N-Telopeptide, nmol BCE/L	5–24	11 (1)	210 (140)	<0.05
Bone AP, U/L	15–41	20 (3)	62 (17)	<0.01
Total AP, U/L	30–130	59 (13)	106 (32)	<0.05
Urea N, mg/L	80–180	160 (20)	630 (110)	<0.001
Creatinine, mg/L	4–12	9 (1)	106 (19)	<0.001
Potassium, mEq/L	3.5–5.0	4.4 (0.4)	5.5 (0.8)	<0.05
Sodium, mEq/L	135–145	140 (1)	137 (3)	NS
Chloride, mEq/L	95–106	106 (3)	100 (4)	NS
CO <sub>2</sub> , mmol/L	24–31	27 (4)	23 (1)	NS
Phosphate, mg/L	25–45	36 (4)	53 (25)	NS
Calcium, mg/L	84–102	93 (2)	91 (4)	NS
Calcium mmol/L	2.1–2.55	2.32 (0.05)	2.27 (0.1)	NS
Albumin, g/L	32–46	41 (1)	35 (2)	<0.01
Glucose, mg/L	700–1100	950 (20)	1140 (200)	NS
Glucose, mmol/L	3.9–6.1	5.3 (0.1)	6.3 (0.1)	NS

<sup>a</sup> NS, not significant; BCE, bone collagen equivalents.

those in a control group. Other reports involving persons who received  $^{41}\text{Ca}$  describe data from single individuals (12, 13, 32). Using a 10-nCi  $^{41}\text{Ca}$  dose, Johnson et al. (32) were able to measure  $^{41}\text{Ca}/\text{Ca}$  ratios for more than 800 days and attributed small fluctuations in  $^{41}\text{Ca}/\text{Ca}$  to menstrual cycle phases. Freeman et al. (12) compared calcium kinetics determined with dual stable isotopes with  $^{41}\text{Ca}$  measurements and showed good agreement. They showed that a surgically menopausal female was losing  $\sim 100$  mg of calcium from bone per day (12). A separate publication by Freeman et al. (13) showed results suggesting that  $^{41}\text{Ca}/\text{Ca}$  measurements had less variability than N-telopeptide and that  $^{41}\text{Ca}/\text{Ca}$  ratios appear to decrease when bisphosphonates are administered (13). Although these earlier studies were limited in size, they

suggest intriguing possibilities for applications of this technology.

It is difficult to prove that the  $^{41}\text{Ca}$  we administered ended up in the bone of the study participants, but several lines of evidence support this conclusion. It is well known that 99% of total body calcium is stored in bone. In addition, animal experiments have demonstrated that 24 h after administration of  $^{45}\text{Ca}$ , an average of 85% of the absorbed calcium is found in the skeleton (38). Using  $^{41}\text{Ca}$ , we have shown that 4 weeks after mice received 0.55 nCi of  $^{41}\text{Ca}$ , the concentration of  $^{41}\text{Ca}$  was 5000 times higher in bone than other soft tissue (39). Bronner et al. (40) showed that 5 days after adolescent boys received  $^{45}\text{Ca}$  intravenously, the total quantity of  $^{45}\text{Ca}$  excreted did not exceed 15% of the dose. These authors contended that nearly all of the calcium that enters the body at a given time is first retained, presumably in the skeleton (40). It should be noted that initially after an intravenous dose of labeled calcium, most of the label is found in quiescent surface that is in dynamic equilibrium with blood (41).

The rapid initial decrease in  $^{41}\text{Ca}/\text{Ca}$  that we observed in plasma of ESRD patients indicates a faster disappearance of  $^{41}\text{Ca}$  from the central compartment in the patients than in the controls. We hypothesize that the increased bone turnover in ESRD patients increases the available sites of mineralization and hence incorporation of the tracer into the bone formation process. This produces a more rapid initial clearance of  $^{41}\text{Ca}$  from the serum compartment. One could also argue that the  $^{41}\text{Ca}$  we administered to the ESRD patients was being cleared during hemodialysis, thus explaining the more rapid decrease of the tracer in these patients. Using  $^{47}\text{Ca}$ , however, Cochran et al. (35) demonstrated that ESRD




Fig. 4.  $^{41}\text{Ca}$  is cleared more rapidly from blood by ESRD patients than by controls.

Kinetics of 10-nCi intravenous  $^{41}\text{Ca}$  dose over 180 days in controls and ESRD patients. ■ and △ represent the mean (SD) values for  $^{41}\text{Ca}/\text{Ca}$  in the controls ( $n = 4$ ) and the ESRD patients on hemodialysis ( $n = 4$ ), respectively. The area under the curve for the ESRD patients was significantly less than that for the controls ( $P < 0.005$ ).

patients with hyperparathyroidism retained 94%–98% of the  $^{47}\text{Ca}$  and that the dialysate and fecal radioactivity were minor components (35). Clearly, more studies are needed with  $^{41}\text{Ca}$ , but this technology offers the potential for monitoring calcium kinetics over a time period of years, using doses of tracer that pose a radiation risk that is too small to measure.

Management of renal osteodystrophy remains a difficult challenge faced by nephrologists. Current treatment guidelines suggest that in ESRD patients, nephrologists should aim to keep the circulating concentration of intact PTH in a range from 2 to 4 times the upper limit of the reference interval, although the form of PTH to be measured and the target concentration range remain controversial (25, 42). The need for improved tools for managing renal patients is clear because the clinical methods currently used to measure renal osteodystrophy are indirect. The ability to directly monitor calcium metabolism may improve our understanding of bone health, which is increasingly recognized as a predictor of long-term outcomes in dialysis patients. In a recent review of the impact of mineral metabolism on mortality and morbidity in >40 000 hemodialysis patients, the attributable risk associated with disorders of mineral metabolism (18%) was greater than that for inefficient dialysis (5%) or anemia (11%) (43).

In summary, we have characterized a dosing material in terms of chemical and radiologic purity in addition to sterility and pyrogenicity. This dose was administered to 8 human volunteers, and no adverse effects were observed. Using the methodology described in this report, we were able to measure  $^{41}\text{Ca}/\text{Ca}$  isotope ratios for more than 170 days in volunteers who received 10 nCi of  $^{41}\text{Ca}$ . The data presented show that this method can be used to monitor differences in calcium kinetics between healthy individuals and patients with ESRD and hyperparathyroidism.

We are grateful to the healthy individuals and patients who volunteered for this study. The trust they showed in taking a experimental radioactive compound, in the hope that they could contribute to improving medical care, reaffirms our belief in the goodness of human spirit. We also thank Kathleen Amarillas and the rest of the staff in the hemodialysis unit at the VA San Diego Healthcare System for their help. Finally, we thank Kunihiko Nishizumi for providing the stock  $^{41}\text{Ca}$  dose material and Todd Wooddy and Mark Stoyer for establishing the radiopurity of the  $^{41}\text{Ca}$ . This work was performed in part under the auspices of the US Department of Energy by University of California, Lawrence Livermore National Laboratory under contract W-7405-Eng-48; it was also supported by the Veterans Administration and NIH.

## References

1. Neer R, Tully G, Schatz P, Hnatowich DJ. Use of stable  $^{48}\text{Ca}$  in the clinical measurement of intestinal calcium absorption. *Calcif Tissue Res* 1978;10:5–11.
2. deGrazia JA, Ivanovick P, Fellows H, Rich C. A double isotope method of measurement of intestinal absorption of calcium in man. *J Lab Clin Med* 1965;66:822–9.
3. Lee WTK. Use of stable calcium isotopes ( $^{42}\text{Ca}$  and  $^{44}\text{Ca}$ ) in evaluation of calcium absorption in Beijing adolescents with low vitamin D status. *Food Nutr Bull* 2002;23:42–7.
4. Vogel JS, Turteltaub KW, Finkel R, Nelson DE. Accelerator mass spectrometry. *Anal Chem* 1995;67:353A–9A.
5. Freeman SPHT, Vogel JS. Biomedical accelerator mass spectrometry. *Int J Mass Spectrom Ion Process* 1995;143:247–56.
6. Vogel JS, Turteltaub KW. Accelerator mass spectrometry in biomedical research. *Nucl Instrum Methods Phys Res B* 1994;92:445–53.
7. Lu CM, Burton DW, Fitzgerald RL, Deftos LJ, Buchholz BA, Vogel JS, et al. Mass spectrometric immunoassay for parathyroid hormone-related protein. *Anal Chem* 2002;74:5507–12.
8. Turteltaub KW, Felton JS, Gledhill BL, Vogel JS, Southon JR, Caffee MW, et al. Accelerator mass spectrometry in biomedical dosimetry: relationship between low-level exposure and covalent binding of heterocyclic amine carcinogens to DNA. *Proc Natl Acad Sci U S A* 1990;87:5288–92.
9. Buchholz BA, Fultz E, Haack KW, Vogel JS, Gilman SD, Gee SJ, et al. HPLC-accelerator MS measurement of atrazine metabolites in human urine after dermal exposure. *Anal Chem* 1999;71:3519–25.
10. Hillegonds DJ, Fitzgerald R, Herold D, Lin Y, Vogel JS. High-throughput measurement of  $^{41}\text{Ca}$  by accelerator mass spectrometry to quantitate small changes in individual human bone turnover rates. *J Assoc Lab Autom* 2004;9:99–102.
11. Paul M, Ahmad I, Kutschera W. Half-life of  $^{41}\text{Ca}$ . *Z Phys A Hadrons Nuclei* 1991;340:249–54.
12. Freeman SPHT, King JC, Vieira NE, Woodhouse LR, Yergy AL. Human calcium metabolism including bone resorption measured with  $^{41}\text{Ca}$  tracer. *Nucl Instrum Methods Phys Res B* 1997;123:266–70.
13. Freeman SPHT, Beck B, Bierman JM, Caffee MW, Heaney RP, Holloway L, et al. The study of skeletal calcium metabolism with  $^{41}\text{Ca}$  and  $^{45}\text{Ca}$ . *Nucl Instrum Methods Phys Res B* 2000;172:930–3.
14. Suarez JA, Rodriguez M, Espartero AG, Pina G. Radiochemical analysis of  $^{41}\text{Ca}$  and  $^{45}\text{Ca}$ . *Appl Radiat Isot* 2000;52:407–13.
15. Demers LM. Biochemical markers in the management of patients with metastatic bone disease. *Clin Chem* 1999;45:1131–2.
16. Kamel S, Fardellone P, Meddah B, Lorget-Gondelmann F, Sebert JL, Brazier M. Response of several markers of bone collagen degradation to calcium supplementation in postmenopausal women with low calcium intake. *Clin Chem* 1998;44:1437–42.
17. Kleerekoper M. Biochemical markers of bone turnover: why theory, research, and clinical practice are still in conflict. *Clin Chem* 2001;47:1347–9.
18. Chailurkit L, Ongphiphadhanakul B, Piaseu N, Saetung S, Rajatanavin R. Biochemical markers of bone turnover and response to bone mineral density to intervention in early postmenopausal women: an experience in a clinical laboratory. *Clin Chem* 2001;47:1038–88.
19. Vesper HW, Demers LM, Eastell R, Garnero P, Kleerekoper M, Robins SP, et al. Assessment and recommendations on factors contributing to preanalytical variability of urinary pyridinoline and deoxypyridinoline. *Clin Chem* 2002;48:220–35.
20. Watts NB. Clinical utility of biochemical markers of bone remodeling. *Clin Chem* 1999;45:1359–68.

21. Trueba D, Sawaya BP, Mawad H, Malluche HH. Bone biopsy: indications, techniques and complications. *Semin Dial* 2003;16:341–5.
22. Nishiizumi K, Caffee MW, DePaolo DJ. Preparation of  $^{41}\text{Ca}$  AMS standards. *Nucl Instrum Methods Phys Research B* 2000;172:399–403.
23. US Pharmacopeia. United States pharmacopeial convention USP 25. NF 20, the national formulary. 85. Bacterial endotoxin test. Rockville, MD: US Pharmacopeia, 2002:1889–93.
24. US Pharmacopeia. United States pharmacopeial convention USP 25. NF 20, the national formulary. 71. Sterility tests. Rockville, MD: US Pharmacopeia, 2002:1878–83.
25. Elder G. Pathophysiology and recent advances in the management of renal osteodystrophy. *J Bone Miner Res* 2002;17:2094–105.
26. Eknoyan G, Levin A, Levin NW. Bone metabolism and disease in chronic kidney disease. *Am J Kidney Dis* 2003;(Suppl 42):1–201.
27. Deftos, LJ. Radioimmunoassay for parathyroid hormone. In: Jaffee W, Behrman H, eds. *Methods in hormone radioimmunoassay*. New York: Academic Press, 1974:231–47.
28. McGonigle RJ, Wallin JD, Husserl F, Deftos LJ, Rice JC, O'Neill WJ Jr, et al. Potential role of parathyroid hormone as an inhibitor of erythropoiesis in the anemia of renal failure. *J Lab Clin Med* 1984;104:1016–26.
29. Lin Y, Hillegonds DJ, Gertz ER, Van Loan MD, Vogel JS. Protocol for assessing bone health in humans by tracing long-lived  $^{41}\text{Ca}$  isotopes in urine, serum, and saliva samples. *Anal Biochem* 2004;332:193–5.
30. Freeman SPHT, Serfass RE, King JC, Southon JR, Fang Y, Woodhouse LR. Biological sample preparation and  $^{41}\text{Ca}$  AMS measurements at LLNL. *Nucl Instrum Methods Phys Res B* 1995;99:557–61.
31. CRC. Table of the isotopes. In: *CRC handbook of chemistry and physics*, 81st ed. Boca Raton, FL: CRC Press, 2000:11–51.
32. Johnson RR, Berkovits D, Boaretto E, Gelbart Z, Ghelberg S, Meirav O, et al. Calcium resorption from bone in a human studied by  $^{41}\text{Ca}$  tracing. *Nucl Instrum Methods Phys Res B* 1994;92:483–8.
33. Cember H. *Introduction to health physics*, 3rd ed. New York: McGraw-Hill, 1996:169–223.
34. National Council on Radiation Protection and Measurements (NCRP). Report no. 94. Exposure of the population of the United States and Canada from natural background radiation. Bethesda, MD: NCRP, 1987:22.
35. Cochran M, Cochran M, Neville A, Marshall EA. Comparison of bone formation rates measured by radiocalcium kinetics and double-tetracycline labeling in maintenance dialysis patients. *Calcif Tissue Int* 1994;54:392–8.
36. Muller P, Bushaw BA, Blaum K, Diel S, Geppert Ch, Nahler A et al.  $^{41}\text{Ca}$  ultratrace determination with isotopic selectivity  $> 10^{12}$  by diode-laser-based RIMS. *Fresenius J Anal Chem* 2001;370:508–12.
37. Moore ID, Bailey K, Greene J, Lu ZT, Muller P, O'Connor TP. Counting individual  $^{41}\text{Ca}$  atoms with a magneto-optical trap. *Phys Rev Lett* 2004;92:153002.
38. Harrison HE, Harrison HC. Studies with radiocalcium: The intestinal absorption of calcium. *J Biol Chem* 1951;188:83–90.
39. Burton DW, Hillegonds DJ, Fitzgerald RL, Herold DA, Deftos LJ. Pharmacokinetics of the  $^{41}\text{Ca}$  distribution in mice [Abstract]. *J Bone Miner Res* 2004;19:S239.
40. Bronner F, Harris RS, Maletskos CJ, Benda CE. Studies in calcium metabolism. The fate of intravenously injected radiocalcium in human beings. *J Clin Invest* 1956;35:78–88.
41. Parfitt AM. Misconceptions (3): calcium leaves bone only by resorption and enters only by formation. *Bone* 2003;33:259–63.
42. Santini SA, Carrozza C, Vulpio C, Capoluongo E, Luciani G, Lulli P, et al. Assessment of parathyroid function in clinical practice: which parathyroid hormone assay is better. *Clin Chem* 2004;50:1247–50.
43. Block GA, Klassen PS, Lazarus JM, Ofsthun N, Lowrie EG, Chertow GM. Mineral metabolism, mortality, and morbidity in maintenance hemodialysis. *J Am Soc Nephrol* 2004;15:2208–18.